

# Flexible Chargebags

The Ezi-Dock Chargebag is a High Performance Flexible Process Package designed to deliver maximum product recovery when handling Pharmaceutical and Chemical Products.



## Product Description

One of the many unique features of the Ezi-Dock Chargebag is the 2 ply body design. The 2 ply design not only adds much needed strength to the Chargebag, when connected to heavy "Passives" but also adds a further layer of protection against the risk of pinholing.

Manufactured from a <USP87> + <USP661> Compliant LDPE with a Permanent Non Migrating Anti-Static Additive.



## Specification

Ezi-Dock Chargebags are available in 6 standard sizes (other sizes are available on request) each can be supplied with either a 2/4/6inch Tri-Clamp connection:-

- 2 Litre
- 5 Litre
- 10 Litre
- 15 Litre
- 25 Litre
- 40 Litre

Chargebags are supplied individually sleeved and packed 25 per carton (52x35x27cm)

## Summary

- <USP87> + <USP661> Compliant LDPE with Permanent Anti-Static Additive
- 2 Ply Body Construction
- 2/4/6/inch Tri-Clamp Ferrule - BS4825-3
- Individual Bag Label Showing Batch & Bag Number
- Available Gamma Sterilised
- Lower Transport & Storage Costs
- HDPE Blanking Caps available
- Handling aids available
- Comprehensive Range of Gaskets/Clamps available



2 inch, 4 inch & 6 inch Tri-Clamp Connections



4inch & 6inch Tri-Clamp Blanking Caps